

ULTRALUX HEALTH




Red Wrap

USER MANUAL

TABLE OF CONTENTS

4	Meet the UltraLux Red Wrap
6-7	Health Benefits
8-9	Package Contents
10	How to Set Up
11	Read Before Use
13-14	Operating Modes
15	Settings & Troubleshooting
16	Optimize Your Results
18	Frequently Asked Questions
19-20	Maintenance & Care
21	Warranty
22	Product Specifications
23	Contact Us



6

Health Benefits

8

Package
Contents



10

How to Set Up



16

Optimize Your
Results



Designed for Targeted Relief, Not Just Surface-Level Comfort

Welcome. You're about to experience a simple, effective way to support your body's recovery and overall wellness.

The **UltraLux Red Wrap** uses advanced red light therapy to deliver targeted support where you need it most — helping your body feel better, faster.

Designed to fit seamlessly into your routine — no complicated setup, no interruption.

Just consistent, effective support, right where your body needs it most.

This guide will walk you through setup, usage, and tips to get the most out of your Red Wrap.



Feel the Difference With

- Targeted pain and inflammation relief
- Accelerated muscle recovery
- Increased blood flow and circulation
- Consistent, daily recovery support
- Supports healthier-looking skin and cellular rejuvenation
- Supports mobility and flexibility
- Supports relaxation and stress relief

HEALTH BENEFITS

The UltraLux Red Wrap supports skin health, muscle recovery, and overall wellness using targeted red light therapy. This guide highlights key benefits and how to get the most out of your sessions.



Increased Collagen

Collagen is a structural protein that provides support to the skin. Red light therapy has been associated with increased collagen production, contributing to the strength and resilience of tissues.



Stimulation of Cellular Activity

Red light therapy enhances cellular function by promoting increased energy production (ATP) within the cells. This may lead to improved cell metabolism and overall cellular health in the body.



Reduction of Inflammation

Red light therapy can reduce inflammation in the skin and muscles. This can be particularly relevant for conditions associated with skin and muscle inflammation in which red light may further support tissue repair and regeneration, aiding in healing.



Increased Blood Circulation

Red light therapy stimulates blood flow in the the body. Improved blood circulation can deliver more oxygen and nutrients to the cells, aiding in the healing process of muscle and also in pain relief.



Mental Health

Exposure to red light is believed to stimulate cellular activity, particularly in the mitochondria, leading to enhanced energy production and improved neurotransmitter function. The release of certain neurotransmitters (such as serotonin and endorphins) as a result of red light exposure can play a key role in mood regulation and emotional well-being.

HEALTH BENEFITS

The UltraLux Red Wrap combines **Red (660nm)** and **Infrared (850nm)** wavelengths which are ideal for providing therapeutic relief.



RED LIGHT 660nm


Treats arthritis and other joint inflammation, aids in the regeneration of stem cells, reduces hypertension, aiding in stroke-prevention, and accelerates muscle building and recovery.



INFRARED LIGHT 850nm

Relief from injury and chronic pain, reduction of subcutaneous fat cells, increases metabolism and digestions, and improves athletic and cardiovascular performance.

NANOMETERS*



Disclaimer: Since red light therapy supports a variety of health conditions, we cannot list every condition, nor can we identify if red light is the solution for your specific condition. If you are curious about how red light might help with your particular health issue, we recommend researching or contacting your doctor. A great place to start is the National Institutes of Health website at nih.gov.




ULTRALUXHEALTH

PACKAGE CONTENTS



1 x UltraLux Red Wrap



1 x USB-C Cable



1 x User Manual



1 x USB-C
Power Adapter


HOW TO SET UP

Before using your Red Wrap, ensure you have:

- All included components
- Access to a power outlet

Step 1: Connect to Power

Plug the provided USB-C cord into the port on the side of the Red Wrap. Then connect the other end to the power adapter and plug it into a wall outlet.




Step 2: Position the Wrap

Place the Red Wrap on the desired area of your body. Ensure the LED side is facing your skin for proper use.




Step 3: Secure the Wrap

Wrap the strap around your body and fasten the velcro to create a snug, comfortable fit. Avoid wrapping too tightly.



Step 4: Start Your Session

Press the power button to turn on the device. Select your preferred light setting (see functions page for details). The Red Wrap will automatically shut off after the session is complete.



READ BEFORE USE

Before using your UltraLux Red Wrap, please read this manual carefully. This device is not intended for use by children, pregnant individuals, or anyone unable to remove the device independently. Consult a healthcare professional before use if you have any medical conditions, are taking medications that cause light sensitivity, or have sensitive skin.

Important Safety Guidelines for Product Use

- This device is not suitable for children or pets
- Use only the provided power supply and accessories
- Do not use the Red Wrap near water or immerse it in any liquid
- Do not use near flammable materials such as aerosols or open flames
- Do not stare directly into the LED light panels
- Ensure all cords and connections are secure before use
- Do not use if the device or cables are damaged
- Limit use to recommended session times (do not exceed 30 minutes)
- Use only on appropriate areas of the body (avoid direct use on eyes or face)

If discomfort occurs during use, discontinue immediately.



OPERATING MODES



All Red Wrap functions are controlled by the **Power Button**  located on the side of the device. The Red Wrap must be **plugged into power** to operate.

When you plug in your Red Wrap, the Power and Settings buttons light up with a soft blue, indicating power but not that the Pulse Therapy Mode is active.

Continuous Therapy Mode

This mode provides uninterrupted light exposure and is ideal for targeting specific areas.

For best results & to optimize your results, follow the instructions on page **16**



To initiate Continuous Mode, follow the steps below:



1x



Plug in the device to a power source & quick-press the Power Button once - the **Power Button will then turn green**. This will activate both Red & Near-Infrared LEDs (660nm & 850nm)

2x



Press again to trigger **only the Red LEDs** (660nm)

3x



To activate **only the Near-Infrared LEDs** (850nm) press the Power Button again. Near-infrared light is invisible to the eye, but the Power Button will remain green to indicate it's active.

4x



Press again to **shut down** all LEDs.

OPERATING MODES



All Red Wrap functions are controlled by the Power Button , located on the side of the device.

The Red Wrap must be plugged into power to operate.

IMPORTANT

When you plug in your Red Wrap, the Power and Settings buttons light up with a soft blue, indicating power but not that the Pulse Therapy Mode is active.



Pulse Therapy Mode

This mode delivers light in pulsed intervals to stimulate blood flow and cellular activity, ideal for **hair growth** and **circulation**.

TO INITIATE PULSE THERAPY MODE:



POWER CONNECTED
Buttons light up green



PULSE MODE ACTIVE
Lights change to blue

After plugging in your Wrap, **quick-press the Power Button until it turns green.**

Then, **press and hold the Power Button for 3–5 seconds** until it changes from green to blue, activating all LEDs (660nm & 850nm) in Pulse Mode.

EASY CONTROL: ONE BUTTON, THREE CLICKS

1x



Red Light Mode

Activates only the **Red LEDs (660nm)**

2x



Near-Infrared Mode

Activates only the **Near-Infrared LEDs (850nm)**
The blue Power Button confirms Pulse Mode is active — even though near-infrared light is not visible.


3x



Power Off

Press again to power off all LEDs.


SETTINGS & TROUBLESHOOTING




Adjust Brightness Settings

There are 2 different Brightness Settings, signified by the **three lines next to the power button**. You can adjust the brightness on both Continuous and Pulsed Therapy Modes.

All Brightness Settings are initiated using the **Settings Button** located next to the power button. The default brightness is 75%.



75% Brightness



100% Brightness



Troubleshooting your Red Wrap

If you are experiencing difficulty with your Red Wrap, we recommend powering off your device and unplugging each cord from the Red Wrap. Then unplug the power cord from the power outlet of your wall.

Then reconnect the cord to the Red Wrap, and the connector to the power outlet. This will support the device in resetting.



***Disclaimer:** When the Infrared is on, you may not see it because near-infrared light is not visible to the human eye. This might give the impression that the LEDs are turned off, however just know they are fully operational.

OPTIMIZE YOUR RESULTS



Inflammation Relief

If you experience symptoms of arthritis and inflammation, or desire to take preventative measures, select Continuous Therapy Mode, activating only the Red light LEDs. Use for one 15-minute sessions every day for a minimum of 4 weeks.

1

SESSION
PER DAY



15 MIN



75%
INTENSITY

4

WEEKS



Weight Loss and Metabolism

For weight loss results, select Pulsed Therapy Mode, activating both the Red light and Infrared light LEDs. Use for one 15-minute sessions every day for a minimum of 4 weeks.

1

SESSION
PER DAY



15 MIN



75%
INTENSITY

4

WEEKS



Pain Relief

If you experience pain relief, select Continuous Therapy Mode, activating only the Infrared light LEDs. Use for one 15-minute sessions every day for a minimum of 4 weeks.

1

SESSION
PER DAY



15 MIN



75%
INTENSITY

4

WEEKS



Muscle Building and Recovery

For muscle improvement and overall recovery, select Pulsed Therapy Mode, activating only the Red light LEDs. Use for one 15-minute sessions every day for a minimum of 4 weeks.

1

SESSION
EVERY DAY



15 MIN



75%
INTENSITY

4

WEEKS



Results may vary. For best results, consistency is key. Always consult with your healthcare professional before starting any new wellness routine, especially if you have a medical condition.



FREQUENTLY ASKED QUESTIONS

How does red light therapy work?

Red light therapy delivers non-invasive light wavelengths that stimulate your cells' mitochondria, boosting energy production and supporting overall cellular function.

Are there side effects from using this device?

Red light therapy is not associated with any adverse side effects if used as directed.

How long will it take to notice results?

While everyone is unique, many users can start to experience pain relief in as soon as one week! We recommend using the Red Wrap daily to experience optimal results.

How long will results last?

If you are consistent with treatment, and maintain an active, healthy lifestyle, your pain relief results could last forever!

Does red light therapy hurt?

No, it should not hurt. Red light therapy should feel like a comfortable, slight warmth. If you experience irritation, stop using and consult a medical professional right away.

What area of the body can be treated?

You can treat any area of the body where you'd like to lose weight or relieve pain. If you have questions about where you can use this red light therapy device, please consult your doctor or medical professional.

Keep Your Red Wrap Performing at Its Best



DAILY CARE

- Wipe clean with a soft microfiber cloth
- Use a gentle non-abrasive cleaner if needed
- Allow the device to fully dry before storing



AVOID

- Do not machine wash
- Do not submerge in water
- Avoid harsh chemical cleaners
- Do not bend or fold tightly around the control unit



STORAGE TIPS

- Store in a cool, dry place
- Disconnect from power after use
- Keep away from excessive heat or moisture



IMPORTANT

Near-infrared light may not always be visible to the human eye, even while the device is functioning properly.

MAINTENANCE & CARE

Keep Your Red Wrap Performing at Its Best



STORAGE:

- Always store flat with the cord detached
- Never fold, roll tightly, or store in a way that creates sharp bends
- Avoid storing coiled or compressed



USAGE:

- When wrapping around a body part, ensure the cord is not under tension or pinched



CORD CARE:

- Never pull, yank, or place weight on the cord
- Avoid bending the cord at sharp angles near the connection points



IMPORTANT TO KNOW:

- LED connection issues are caused by bending, folding, or cord tension over time
- Once LEDs go out in an area the damage is not repairable
- This is the most common reason for wrap failures and is not covered under warranty if caused by improper storage or use



Proper care today ensures your Red Wrap continues to deliver consistent, effective performance.

WARRANTY

Your UltraLux Red Wrap includes a 1-year manufacturer warranty for the original purchaser. To qualify, the Red Wrap must be in like-new condition aside from the reported defect.

To start a warranty claim, contact support@ultraluxhealth.com with photo or video evidence of the issue. UltraLux Health will review the claim and determine coverage at its discretion. If approved, your Red Wrap will be repaired or replaced at no cost. Valid warranty claims cover labor, defective parts, and return shipping of the repaired or replacement product during the first year of coverage.

For video instructions, scan this QR code
with your phone's camera



PRODUCT SPECIFICATIONS

Product Name	UltraLux Red Wrap
Rating Power	38 Watts
Actual Power	15 Watts
Product Size	45x0.19x5.25in
LED Quantity	75 diodes / 150 individual LEDs (2 LEDs per diode)
Net Weight	0.29 lb
Wavelength	660 nanometers ; 850 nanometers
Irradiance	0 inch: > 100 mW/cm ²
Timer	20 minutes
Material	Fabric
Lifespan	50,000 hours

CONTACT US

Email

support@ultraluxhealth.com

Website

ultraluxhealth.com

Phone

385-360-9755

Hours

Monday – Friday

8:00 AM – 4:00 PM MST